

The versatile power wheelchair.

CentroGlide is our leading model. It's a robust mid-wheel-drive power wheelchair, designed and built by us in Australia to provide remarkable stability and trouble-free service for years.

Completely customisable.

CentroGlide is a fully-scripted wheelchair, built to the occupant's precise requirements.

Independence and freedom of movement.

CentroGlide has been designed to allow remarkable manoeuvrability in tight spaces, and stability over kerbs. It comes with 50° seat tilt as standard, with options for vertical lift and anterior reach to allow the users to reach high shelves or transfer from the chair in a standing position.

Power through the day.

CentroGlide is built with the functionality and reliability to give users total mobility, indoors or out.

Designed for comfort and ease of use.

We have put decades of user feedback and experience to develop all the features that make CentroGlide a popular choice.



Customised seating.

For ultimate comfort many seating options are available. Choose customised cushion, backrest and headrest.



Active suspension.

Our patented 6-wheel APT suspension improves traction and aids ascending and descending kerbs safely without affecting the stability and angle of the seat.



Wheels.

CentroGlide has 14" drive wheels, 8" and 7" castor wheels as standard. Other sizes are available.



Colour choices.

Choose from an extensive colour palette for the ultimate customisation.



Electronics.

Fully programmable R-net joystick control system. Compatible with many other 3rd party specialty electronics.



Power seating.

As standard, the seat has a 50° tilt. Options include power seat lift and anterior reach.



Battery range.

The powerful battery gives you a range of up to 40km on a single charge.



CentroGlide sets the benchmark.

Features.

- Power base with removable seating module
- High quality German 2 pole with optional 4 pole high torque motors
- Vehicle tie down points
- Low maintenance with accessible electronics
- Turning circle between walls 1145mm
- Active Pozi-Track (APT) kerb climbing ability
- 80mm ground clearance
- Max ramp incline/decline 10° (as per AS3695 testing)
- 4.1/3.50-6 (14") drive wheels
- Electronic brakes with manual park
- R-net LCD Std controls

Safety.

- Full compliance AS 3695 175kg
- Crash tested to AS 3696:19 to max user weight 175kg

Arrange a test drive.

Chat to a Glide Wheelchair specialist and experience the extraordinary features of the GlideCentro in person. Find your nearest distributor at glide.com.au/nearest

Specifications.



Recommended User Weight
Up to 175kg AS 3695 full compliance



Overall Length
900mm (foot plate removed)



Seat Widths
31 to 61cm



Overall Width
620mm (wheel base)



Seat Depths
31 to 58cm



Seat Height
46cm (top of rail using 0-46 set up)



Back Height
1030mm Top of push handle (custom height available)



Arm Rests
Flip up (adjustable)



Battery
Up to 40km range



Leg Rests
Swing away detachable



Frame Colour
Wide choice



Brakes
Electronic and manual park



www.glide.com.au

(08) 9345 3400 | sales2@glide.com.au

25 Ledger Road, Balcatta, Western Australia 6021

@GlideProductsAustralia Glide Products